



Disposable Diagnostic Pen Torch With Pupil Gauge



6 Pens Per pack
High Intensity

www.med-supply.com.au Reorder Code: DT6205



Disposable Penlight Torch

Self focusing, high intensity bright light without the inconvenience of changing batteries and bulbs. Ideal for the physician, nurse and EMT and can be easily clipped to pockets. All penlights include two "AAA" batteries as a sealed unit and high intensity pre-focused lens and bulb.

Product Product Code Packaging

Disposable Diagnostic Penlight Torch	DT6205	6/Bag, 36/Box
--------------------------------------	--------	---------------

- Completely disposable, single patient use
- Suitable for testing pupillary dilation
- Includes pupil gauge scale measuring from 2mm - 9mm
- Single click clip activated on/off switch
- Concentric beam bulb with pocket clip attachment
- Sealed unit with 2 x AAA size batteries included
- Tamper proof design
- Sold in pack of 6