

# Reusable Antibacterial Mask

## 3 Layers



### Outer Layer

Water & splash resistant  
UPF 50+ sun protection  
100% polyester

### Middle Layer

3 filter function (air & dust)  
100% polypropylene

### Middle Layer

99.9% antibacterial nanofiber  
65% polyester - 35% cotton



Liquid  
Repellent



Kills 99%  
Bacteria



Reduces  
Odour



30X  
Washable



UPF  
50+  
UPF 50+  
Rated

## Instructions

**Important:** Wash hands before and after handling the mask. Apply the mask using the two elastic bands. Place mask in front of nose and mouth and attach the elastic bands behind the ears. Avoid touching and sliding the mask. Do not place the mask on forehead or chin. Remove the mask using the two elastic bands and store or wash as needed.

## Features

Water and splash resistant fabric  
99.9% antibacterial protection\*  
UPF 50+ sun protection\*  
Individual sterile packaging  
Adapts to any face shape

\*As tested and certified by Intertek

## Washing

Quality tested for 30 washes.  
Wash at 60° and air dry.  
Do not iron on high heat.  
Dispose after 30 washes

Code: L0Y4B

Size: Adult

Colour: Black

Code: L0Y4W

Size: Adult

Colour: White

Code: L0Y4DB

Size: Adult

Colour: Blue

## DISTRIBUTOR:

Medsupply Australia Pty Ltd  
153 Victor Road, Sydney 2099  
TEL: 1300 587 200  
www.med-supply.com.au



STERILEEO



1 234567 890128 >