



# COVID-19 Antigen Rapid Test Cassette




## Nasal Swab

### (For Self-Testing)

The COVID-19 Antigen Rapid Test Cassette (Nasal Swab) is a rapid chromatographic immunoassay for the qualitative detection of COVID-19 antigen in Nasal Swab. The identification is based on the monoclonal antibodies specific for the Nucleocapsid (N) protein of SARS-CoV-2. It is intended to aid in the rapid differential diagnosis of current COVID-19 infections.

#### **Painless**

 Avoid the discomfort of nasopharyngeal swab sampling

#### **Convenient**

 Easy to operate, Interpret result at 10 minutes

#### **Accurate**


 Excellent performance compared to molecular methods



For customer service, call 1300 587 200  
For all order enquiries, email [orders@med-supply.com.au](mailto:orders@med-supply.com.au)

©2012 Medsupply Australia Pty Ltd. All rights reserved. A.B.N. 33 158 708 975






## ► MATERIALS

- Test cassette
- Prefilled Buffer
- Dropper tip
- Sterile Nasal Swab
- Package Insert
- Quick Reference Guide
- Disposal bag

## ► Convenient Procedure



## ► Ordering Information

Cat.No.	Product	Specimen	Format	Pack
ICOVN-C82H	COVID-19 Antigen Rapid Test Cassette (For Self-testing) CE1434	Nasal Swab	Cassette	2T
ICOVN-C85H	COVID-19 Antigen Rapid Test Cassette (For Self-testing) CE1434	Nasal Swab	Cassette	5T
ICOVN-C810H	COVID-19 Antigen Rapid Test Cassette (For Self-testing) CE1434	Nasal Swab	Cassette	10T
ICOVN-C820H	COVID-19 Antigen Rapid Test Cassette (For Self-testing) CE1434	Nasal Swab	Cassette	20T